INNev

### INIR30



10

INIR36



30 nos. of dropping bottle + basket for caps

36 nos. of 250ml erlenmeyer flasks, volumetric flasks, beakers, test tube, etc. (distance between injectors 99\*90mm)

INIR38




38 nos. of cylinder. To simultaneously wash mearsure cylinder 16 nos. of 5~25ml, 12 nos. of 50~250ml, 6 nos. of 500ml and 4 nos. of 1000~2000ml

INIR60



60 nos. of erlenmeyer flasks, volumetric flasks, beakers, test tube, etc. (distance between injectors 72\*70mm)

**INIR80** 



80 nos. of small beaker rack. To wash small beakers, centrifuge tubes, volumetric flasks, etc. (distance 56mm)

### **INIR11218**



112 nos. of injector for beakers, erlenmeyer flasks, volumetric flasks etc.(distance between injectors 30\*30mm)+18 nos. of injector for beakers, erlenmeyer flasks, volumetric flasks etc. (distance between injectors 90\*86mm)

## **INIR120**

INIR222





INIR260



INAB40



To wash 40 nos. of different specification of absorbin bottle and sample bottle



120 nos. injectors. To wash 120 nos. of test tubes, 120 nos. of less than 50ml volumetric flasks, etc. (distance between injectors 53\*48.5mm)

222 nos. injectors, 222 nos. of autosamper vials, test tubes (distance between injectors 31\*23mm)+1 set of basic basket (bottle caps, stoppers etc.)

260 nos. of autosamper vials, test tubes (distance between injectors 38\*23mm)

# INNZVA

# INRTR20



12

To wash 20 nos. of rubber tube with its inner diameters 3~8mm and length less than 3m

INPR60



60 pieces of pipette. To wash 60 nos. of 550mm pipettes/straws. It can also wash some wide mouth glassware such as beaker

INPR6060



60 nos. of pipette and 60 nos. of autosamper vials

INPDR9070U



70 nos. of culture dishes within diameters 90mm (upper layer)

INPDR9070L



70 nos. of culture dishes within diameters 90mm (lower layer)

#### INMR04/INMR05



INBR01







Medical device rack. 4 or 5 layers, 2 baskets per layer

Without injectors, to wash wide mouth flasks such as beaker, the area is 596mmx513mm.